



System requirements for ProgRes® USB 2.0 microscope cameras

Requirements for the installation of **ProgRes CapturePro software** (included in delivery of ProgRes microscope cameras).

Beside of the latest software version on the CD / memory card, included in delivery of the cameras, we offer a free download after registration.

Download: [ProgRes CapturePro - Software](#)

RECOMMENDED SYSTEM REQUIREMENTS

Personal Computer	Intel i7 (Quad-Core) Processor or comparable 8 GB RAM
Operating System	Win 7 SP1 / Win 8.1 / Win 10 – 64 Bit
Data Interface	USB 3.0 integrated USB 2.0 integrated
Monitor Resolutions	1920x1080 or higher
Graphic Interface	Graphic card equipped with on-board memory

MINIMUM SYSTEM REQUIREMENTS

Personal Computer	Intel i3 (Dual-Core) Processor or comparable 4 GB RAM
Operating System	Win 7 SP1/ Win 8.1 / Win 10 – 32 & 64 Bit
Data Interface	USB 3.0 over PCI Express V1.1 USB 2.0 integrated
Monitor Resolutions	1280x720 or higher